



YURATEK

DURU-II

A professional Approach to Refraction Units

www.yuratek.com

DURU-II

Refraction Units

TECHNICAL SPECIFICATIONS UNIT AND CHAIR

•Unit	DURU-II
•Power	110 or 240 VAc 50-60 Hz
•Electrical Classification	Class I ,Type B
•Table Top Electrical Output	110 or 220 VAc & 0-6/12(SL)
•Head Lamp Dimmable	Led Lamp
•Maximum Load of Table	20 kg + 20 kg
•Chair Carrying Capacity	150 kg
•Chair Stroke	180 mm

STANDARD FEATURES

- Sliding Table for 2 Instruments
- Easy change to Left or Right Version at site
- Electromagnetic Brake in Every Angle
- Chart Projector Plate
- Safety Knee Switch
- Dimmable Slit Lamp Control
- Standard Trial Lens Drawer
- Backlight Illuminated Keypad
- Energy Saving Mode
- Mini Keypad on Table
- Reading Lamp (Near Vision)
- Column with Adjustable LED Lamp
- Height-adjustable Projector Support
- Foldable Armrests
- Chair Rotation : 360°

CONTROL PANEL


- Power On/Off Control
- Chair Up/Down Elevation Control
- Aux On/Off Control
- Projector On/Off Control
- Head Lamp Dimmable Control
- Slit Lamp Illumination Light Control
- Elevated Table Control
- Motorized Phoropter Arm Control
- Electro Magnetic Brake Control
- Reading Lamp (Near Vision) Control

OPTIONAL ACCESSORIES

- Manuel Phoropter Arm
- Autobalance Phoropter Arm
- Sliding Manuel Phoropter Arm
- Motorized Phoropter Arm
- Common Chinrest
- Foot Controller Pedal For Magnet
- Foot Controller Pedal For Chair
- Monitor Support Arm
- Mobile Triple Trail Lens Drawer
- Elevated Table Top
- Double Column
- Wheelchair Access Type
- Inclinable Chair Backrest 45°
- Inclinable Chair Backrest 180°
- Manual Horizontal Chair Movement
- Ophtalmoscope Retinoscope Holder
- Chair Stroke 200/300/400 mm
- NEW** Fixed Triple Trail Lens Drawer
- NEW** Optional Touch Screen
- NEW** Maximum Load of Table(Total 75 kg)
- NEW** Motorized Horizontal Chair Movement (Only with Touch Screen)

Please contact us for customised options, accessories and detailed information.


STANDARD TYPE TECHNICAL DRAWING



Wheelchair Access Type



Double Column Type



Triple Trail Lens Drawers Type



Left or right sided table installation feature allows you to easy adjust the unit's table during the installation.



Backlight Keypad



Optional Touch Screen

YURATEK TIBBİ CİHAZLAR LTD. ŞTİ.

Saray Mah. 47.Sok. No:1 Kahramankazan - ANKARA /TÜRKİYE

Tel : +90.312 395 29 96 Fax : +90.312 395 29 97 e-mail: info@yuratek.com www.yuratek.com

YURATEK



IMPORTANT Yuratek Ltd. Şti. subject to change in design and / or specifications without advanced notice. For detailed info please contact your dealer.